

THE PERFECT PAIR



EZ11Plus

The Perfect Primary Unit

Designed to meet the most current sterilization standards ANSI/AAMI ST55

Features:

- Closed-door HEPA filtered air-drying system to maintain sterility and ensure efficient drying of packs and pouches, this ensures that the door remains closed throughout the dry cycle.
- 11" chamber and 19.8" deep, 5 extra-large trays for pouches and packs
- 4 FDA cleared sterilization programs. Documenting and storing cycle data can be transferred to a PC from the USB port.
- 2 year parts and labor warranty, 10-year chamber warranty

Small, Fast, Economically Priced



Valueklave™ 1730

The Perfect Cassette Alternative Unit

Features:

- Fast Cycle Time, Sterilizes unwrapped instruments in 11 min.*
- Bigger chamber size than a single cassette unit.
- 3 small stainless steel trays no need for special cassettes
- Effectively dries wrapped instruments
- Overall dimensions 17.9"D X 17.4"W x 12"H

*Come up time, sterilization exposure and exhaust

Buy a New EZ11Plus and Valueklave 1730 combo from October 1 thru December 31, 2017 and receive a \$600 manufacturer's rebate

To redeem, the end user must visit www.tuttnauerUSA.com/Q4-2017-combo-rebate, upload the dealer invoice dated October 1 thru December 31, 2017 (with doctor's and/or practice/facility name and date purchased), complete the information required on the form (which includes the serial numbers of both autoclaves), and select "Submit". All information will be verified by Tuttnauer USA upon submission. A confirmation email will be sent for your records. Keep the confirmation until rebate is received. Redemption deadline is January 31, 2018. Please allow 4 weeks for receipt of rebate. Excludes reconditioned and scratch and dent machines. BOTH MACHINES (EZ11PLUS AND VALUEKLAVE 1730) MUST BE PURCHASED AT THE SAME TIME AND BE ON THE SAME INVOICE TO BE VALID. CANNOT BE COMBINED WITH ANY OTHER MANUFACTURER OFFERS OR PROMOTIONS



Tel: (800) 624 5836, (631) 737 4850
Email: info@tuttnauerUSA.com,
www.tuttnauerUSA.com

